Train-On-Request: an On-Device Continual Learning Workflow for Adaptive Real-World Brain-Machine Interfaces


Lan Mei¹, Cristian Cioflan¹, Thorir Mar Ingolfsson¹, Victor Kartsch¹, Andrea Cossettini¹, Xiaying Wang¹, Luca Benini^{1,2}

¹Dept. ITET, ETH Zürich, Switzerland ²DEI, University of Bologna, Italy

BiOCAS 2024 § Biomedical Circuits and Systems Conference October 24-26, 2024 | XI'an China

Poster No. B2P-C01





Continual On-Device Learning (ODL) with short calibration time is needed! Results More successful sessions Unsuccessful / Successful TOR Session observed in last 3 sessions Finetune Finetune Acc.: 93.3% Acc.: 56% (%) (%) Acq. time: Acq. time: 8.3min .6min Adapt. time: Adapt. time: 2 3 4 5 6 7 8 9 10 5 6 7 8 9 10 16.2s 3.2s Subsession Subsession TOR Performance over Multiple Sessions Optimal chain-TL baseline: 60/40 train/test, 86.2% acc., ~10min acq. time per session TOR-TL TOR-LwF TOR-ER #Train. Trials & Acq. Time: TOR methods -> fewer #training trials, less acquisition time TOR-ER -> best among all, #trials 46% fewer than chain-TL, ■-- Chain-TL —— TOR-TL —— TOR-LwF —— TOR-ER acq. time 5.3min per session Accuracy:

- TOR-ER -> best among TOR
- TOR accuracy increases over time, TOR-ER comparable with chain-TL, especially for the last 3 sessions (88.8%)
- 7.6 bit/min ITR for TOR-ER

Accuracy threshold T_{Acc}:

 T_{Acc} increases -> Required #training data increases, accuracy increases in most sessions


#Trials per subsession trls:

 trls = 10 -> an optimal accuracy-training time trade-off, better than trls = 5 and trls = 20

~1mJ, 21.6ms per training step

Adaptation time: 3.2s for 15 epochs << 1.6min acquisition time per subsession

On-site adaptation for an improved user experience!



Conclusion

- Novel TOR ODL workflow introduced for EEG-based systems.
- Tackling inter-session variability through near real-time data and **online finetuning** for out-of-distribution acquisition sessions.
- Achieved a 46% reduction in calibration time, lowering it to 1.6 minutes, with up to 92% accuracy.
- Deployed on GAP9 MCU with 50 mW power, 1 mJ per learning step, and 21.6 ms runtime.
 - A robust, user-friendly, and energy-efficient solution to increase user acceptance of BMIs in non-clinical settings.

Acknowledgment

This work was supported by the Swiss National Science Foundation under grant n. 193813 (PEDESITE) and grant n. 207913 (Tiny-Trainer), by the ETH-Domain Joint Initiative program (project UrbanTwin), and by the ETH Future Computing Laboratory (EFCL).



 T_{Acc} | 70% | 80% | 90% | 100%

trls 5 10 20

Session

On-device Measurements

800mV, 370MHz on GAP9