

Finetuning and Quantization of EEG-Based Foundational Biosignal Models on ECG and PPG Data for Blood Pressure Estimation

Integrated Systems Laboratory (ETH Zürich)


Thorir Mar Ingolfsson, Bálint Tóth, Dominik Senti, Jeffrey Zweidler, Alexandre Elsig, Luca Benini, Yawei Li

thoriri@iis.ee.ethz.ch

COPENHAGEN




International Conference of the IEEE Engineering in Medicine and Biology Society




The Wearable Revolution & Its Processing Gap



Rapid Growth in Wearable Health Monitoring



Heart Rate (ECG,PPG), Brain Signals (EEG) etc.



Wearables gather data and send away






Limited On-device processing!





A Widening Chasm: Powerful AI vs Tiny Devices







How do we push these large brains onto small battery powered edge devices?





Algorithms for HW and HW for Algorithms







FC Domain

Cluster Domain

L2
1.5
MB

1x FPU Flash
2 MB

Cluster Domain

Cluster Domain

Ax FPU

L1 128 kB

9x Core

4x FPU

GAP9





Key limitations:

Memory: L1 128kB

L2 1.5MB

Speed: 400MHz

9 Core

computation

cluster











The Next Frontier: From Specific to Universal Models





X: Raw Data Cheap

Heart Rate (ECG)

~ 10TBs (20k+ hours).

Brain Signals (EEG)


Muscle Activations (EMG)


Y: Labels

Atrial Fibrillation, Seizure, etc.

Expensive

Supervised Learning





Foundation






How do we train on unlabeled data?



LLMs are often pretrained in an intuitive way.

The ocean was very calm and the waves crashed into the rocks





Other methods: Contrastive Learning



CEReBrO's Alternating Attention to the Rescue









Intra-Channel Attention

Self-attention between tokens of the same channel




Inter-Channel Attention

Self-attention between tokens from the same timestep



Biosignal Cross-Modality Knowledge Transfer







Cross-Modality knowledge Transfer brings SoA Results



VitalDB MIMIC-III

Method	Diastolic Blood Pressure MAE	Systolic Blood Pressure MAE
ResUNet+Attention	2.95	4.64
CiGNN	2.79	4.15
Fusion+Attention	3.76	5.32
Our CEReBrO	1.92	3.14

Method	Diastolic Blood Pressure MAE	Systolic Blood Pressure MAE
ResUNet+Attention	1.13	4.55
ResNet+LSTM	2.23	4.12
Conv-LSTM	3.29	4.30
Our CEReBrO	1.57	2.72

Lower → better


SoA results → Cross Modality Knowledge Transfer works!


Best CEReBrO models quantized down to around 3.8 MB → Fits on edge devices




Bringing Foundation Models to the Edge











Specific kernels written to handle alternating attention







Cross-Modality and State-of-the-art Results


First EEG-based Foundation Model used for ECG/PPG usage

→ Cross-Modality Knowledge Transfer



State-of-the-art accuracy diastolic BP (MAE 1.57 mmHg) and surpasses by 1.5× the accuracy of prior works for systolic BP (MAE 2.72 mmHg).

Contact info





Thank you!

